GRiPKIT

BY WEISS ROBOTICS

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TM

1 Introduction

GRIPKIT CR PLUS is Weiss Robotics' intelligent gripping solution for robot arms of Techman Robot. GRIPKIT CR PLUS is available with different gripping modules and is mechanically and electrically fully compatible of Techman Robot.

This operating manual describes the installation, initial setup and operation of the GRIPKIT CR PLUS. It does not replace the operating manual of the used gripper which can be found on the included USB flash drive or Online at www.weiss-robotics.com/gripkit-cr-plus.

1.1 Related documents

Related documents for operating the gripping module, technical drawings and 3D models are available at www.weiss-robotics.com.

Please refer to our General Terms and Conditions for additional information on the warranty, available at www.weiss-robotics.com/gtc.

1.2 Target audience

This manual is intended for integrators and operators of robot applications and all persons maintaining and working with this product. It is advised to keep this and related documents available to the staff at all times. Please read this document carefully and observe in particular the safety and warning notices in this document.

1.3 Notation and symbols

For a better understanding, the following symbols are used in this manual.

Functional or safety-relevant information. Non-compliance may endanger the safety of personnel and the system, damage the device or impair its function.

Additional information for a better understanding of the described facts.

Reference for further information.

2 Basic safety instructions

2.1 Intended use

GRIPKIT CR PLUS was developed for reliably gripping and holding of work pieces or other objects and is intended for mounting on a robot arm. The requirements of the applicable directives as well as the assembly and operation instructions in this document must be observed and adhered to. GRIPKIT CR PLUS may only be used within its defined operating limits and only for industrial applications.

Any other use is considered to be improper, e.g. if the gripping module is used for pressing, curving, lifting, punching or as a clamping or guidance aid for other tools. The manufacturer is not liable for any damage resulting from this improper use.

2.2 Environmental and operating conditions

GRIPKIT CR PLUS and the included gripping module may only be used in the context of its defined application parameters. Make sure, that the module and the gripper fingers are at sufficient size for the application. The environment has to be clean and the ambient temperature must not exceed the specifications in this document. In addition, the environment must be free from splash water and vapors as well as from abrasion or processing dust, except the used gripping module is especially designed for contaminated environments.

2.3 Controlled production

GRIPKIT CR PLUS represents the state of the art and the recognized safety rules at the time of delivery. However, it can present risks if, for example:

- GRIPKIT CR PLUS is not used in accordance with its intended purpose
- GRIPKIT CR PLUS is modified, or not installed or maintained properly
- The EC Machinery Directive, the VDE directives, the safety and accident-prevention regulations valid at the usage site, or the safety and installation notes are not observed

2.3.1 Protective equipment

Provide protective equipment per EC Machinery Directive where required.

2.3.2 Constructional changes, attachments or modifications

Additional drill holes, threads or attachments that are not offered as accessories by

Weiss Robotics may be attached only with written permission of Weiss Robotics.

2.3.3 Special standards

The following harmonized standards are adhered to:

- RF interference, noise and emission according to IEC/EN 61000-6-3
- Fast transients (bursts) on supply and I/O lines according to IEC/EN 61000-4-4
- RF power input on supply and I/O lines according to IEC/EN 61000-4-6
- RF irradiation according to IEC/EN 61000-4-3
- Interference emission according to IEC/EN 61000-6-4 Class A
- Magnetic field with energy frequency according to IEC/EN 61000-4-8
- Electrostatic discharge according to IEC/EN 61000-4-2
- IO-Link communication standard according to IEC 61131-9

2.4 Personnel qualification

The assembly, initial commissioning, maintenance, and repair of the GRIPKIT CR PLUS may be performed only by trained specialist personnel.

Every person called upon by the operator to work on the GRIPKIT CR PLUS must have read and understood the complete operating manual. This applies particularly to occasional personnel such as maintenance personnel.

2.5 Safety-conscious working

Observe the safety and accident-prevention regulations valid at the usage site.

Do not reach into the open mechanism or the movement area of the gripping module. De-energize and set the gripping module pressure-free before performing any work on it.

3 Warranty

The warranty is 12 month from the ex-works delivery date in accordance with the intended use in one shift operation or 20 million gripping cycles and observing the prescribed maintenance and lubrication intervals. Maintenance and lubrication intervals can be found in the operation manual of the included gripping module. Workpiece contacting parts and wearing parts are not covered by this warranty. Refer to our Terms and Conditions for details. GRIPKIT CR PLUS is considered to be defective, if its basic function gripping is no longer given.


4 Model variants

GRIPKIT CR PLUS is available in different model variants.

GRIPKIT CR PLUS	S	L
Gripping module (Part No.)	CRG 30-050 (5010016)	CRG 200-085 (5010017)
Active mounting plate (Part No.)	SF-ISO50-CRG30 (5020061)	SF-ISO50-CRG200 (5020062)

5 Assembly

The gripping module is mounted to the robot arm using the active mounting plate contained in the GRIPKIT CR PLUS. The necessary mounting material (screws, pins, etc.) is included in the GRIPKIT CR PLUS.

Only work on the electrical system of the robot when it is switched off.


Do not kink the cable at the connector! Damage to the connector possible!

With GRIPKIT CR PLUS S: Install the enclosed spacer plates between the flange and the robot. Use the long screws with a diameter of 6 mm for this.

6 Technical specification

PIN	FUNCTION
1	IN2
2	IN1
3	INO
4	Reserved
5	Reserved
6	GND
7	+24V/1.5A IN
8	OUT2
9	OUT1
10	OUT0

Power Consumption

Model	Holding at 25%	Holding at 100%
GRIPKIT CR EASY	0,05A	0,25A

6.1 Input Circuit

The inputs of the gripper module support both NPN and PNP signals. The passive state is detected when the input is open (inactive). The active state is detected when the input is actively pulled to either 0V (NPN) or 24V (PNP).

Input	Function GRIPKIT CR PLUS
INO	GRIP active: grab workpiece passive: workpiece released
IN1	INDEX.0 selects the grip to execute
IN2	INDEX.1 selects the grip to execute

Grip	IN1	IN2
0	p	p
1	a	p
2	a	p
3	a	a

6.2 Output circuit

The gripping state is output via the outputs of the gripping module. The outputs are designed as a push-pull circuit

Out 2	OUT1	OUT0	Gripping state
0V	0V	0V	wait for command (IDLE)
0V	0V	24V	release workpiece (RELEASED)
0V	24V	0V	no workpiece detected (NO PART)
0V	24V	24V	holding workpiece (HOLDING)
24V	0 oder 24V	0 oder 24V	Error

7 Scope of delivery

The following components are included in the GRIPKIT CR PLUS package for UR:

- Gripping module (see chapter „Model variants“)
- Active mounting plate (see chapter „Model variants“)
- Connection kit GRIPKIT CR EASY/PLUS (Part No. 5070018)
- For GRIPKIT CR PLUS S: two spacer plates (Part No. 5020063)
- USB flash drive with software (Part No. 5090008)
- Quickstart GRIPKIT CR PLUS (Part No. 5080030)

Additional Accessories

GRIPKIT CR PLUS does not contain any fingers, since those are usually highly application-dependent.

Order matching fingers from our standard product range separately.

7.1 Content of the included USB flash drive

The included USB flash drive contains a software as well as sample programs and the documentation of the gripping module.

8 Technical data

GRIPKIT CR PLUS	S	L
Gripper type	parallel	
Gripping force	15 - 30 N	75 - 200 N
Workpiece weight (form fit/force fit)	2.7 kg 0.15 kg	4.3 kg 0.15 kg
Full stroke	50 mm	85 mm
Current consumption	350 mA	850 mA
Weight	0.46 kg	1.3 kg
Parameterization	2 grips (position, force)	

Exceeding the specified technical data may reduce the life of the gripping module or permanently damage it. If in doubt, please clarify your application with your responsible sales partner.

9 System requirements

GRIPKIT CR PLUS is compatible to the cobots of Techman Robot. To mount the GRIPKIT CR PLUS size S , an additional spacer pate (Part no. 5020063) is required. The spacer plate is not included in the delivery.


Order matching spacer plates from our standard product range separately.

The following software requirements apply to the operation:


- TM flow: System version 1.72.3550 or newer

10 Device Configurator


To configure your gripper module, you can alternatively use the "Device Configurator" program. You will find this program on the enclosed USB flash drive. Do the installation on a desktop device and open the program. Connect the Gripper to your computer with the enclosed USB cable.

Start the "Device Configurator" program and connect the gripper via "Connect".

Next, activate the control of the gripper and carry out a reference run. Press "Activate" and then "Reference".


Now the gripper is ready for use. You can test the adjusted grips 0-3 by pressing "Grip" and "Release". The respective grip parameters are defined with "Limit value no part" and "Limit value part release".

10.1 Grip direction

The gripping direction is determined by the relation of the position limit values to one another. If the „No Part Limit“ is larger than the „Release Limit“, the gripping module performs an internal grip, if the „Release Limit“ is larger than the „No Part Limit“, an external grip will be performed.

External and internal gripping example

10.1.1 Gripping force

GRIPKIT CR PLUS has an integrated gripping force control. The gripping force can be specified as a percentage of the nominal gripping force of the gripping module. Possible values are 25 to 100 % of the nominal gripping force.

10.2 Define a grip

To define a grip, go to the "Gripping Settings" tab. Select the parameters for "Limit value no part" and "Limit value part release" and confirm the entries with "Apply".

When gripping and releasing, electrical gripping modules stop at the „No Part Limit“ and the „Release Limit“.

You must then specify the gripping force as a percentage of the nominal force. The nominal force of the gripping module used can be found in chapter 8 Technical data. With the parameters "Override gripping" and "Override enable" you can readjust the gripping speed. Both values are set to 100 % by default. When gripping, this value refers to the optimal gripping speed calculated by the gripper depending on the part size and gripping force. You can set values from 1-200 % here. When releasing, 100 % is the absolute maximum speed of the gripper. Values of 1-100 % are possible here.

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
10.3 Set up light ring

The gripping modules have a glare-free LED light ring that can be used to visualize various states of the application. Eight visualization patterns can be parameterized for this purpose. The animation, colour and speed can be changed. To do this, open the "LED" tab and select the various gripping states via "Animation". Confirm the input with "Activate".

10.4 Diagnostics


Under the "Diagnostics" tab you will find the current data on the status of the connected gripping module.

11 Function block installation


The control unit of the GRIPKIT CR PLUS is installed by installing the function blocks on the robot controller. For the installation of the function blocks, an additional NTFS-formatted USB flash drive with the name "TMROBOT" is required.

11.1 Preparation of the USB flash drive

Since the directory structure of the USB flash drive depends on the individual serial number of the robot, it makes sense to export an existing project from the robot to the USB flash drive first. The required directory structure is created automatically. For more information about exporting projects, see the documentation of the robot system. The USB flash drive should now have the following directory hierarchy:


The name of the middle directory depends on the serial number of the robot. Now create a new directory "ComponentObject" and copy the three function blocks "ChangeGrip", "Grip" and "Release", which can be found on the provided USB flash drive, into this directory.



Then remove the USB flash drive and plug it into a free USB port of the robot system.




11.2 Installation on the robot

Select the menu item "System" in the user interface of the robot system, to access the system settings.



Select "Import / Export" in the system settings to import the function blocks from the USB flash drive to the robot system. To import the function blocks, follow the steps below:

At the end of commissioning, the imported function blocks must be activated. Under the menu item "Settings", open the robot settings and activate the function blocks under the "Component" tab by clicking on the "Enable" field. Then save the component list.





12 Programming

After the installation you will find the three new function blocks "ChangeGrip", "Grip" and "Release" in the overview under "Project".


12.1 Grip selection (ChangeGrip)

Regardless of the size, each gripping module supports up to four different grips. With the "ChangeGrip" block one of these four pre-parameterised grips can be selected. To do this, drag the function block to the desired position in the programming field using drag & drop. Then follow the steps below to configure the function block:

Use the "int / ChangeGrip1_var_Preset =" field to set the number of the desired grip (between zero and three) and confirm with "OK".


ENGLISH

TM

12.2 Grip a part (Grip)


Drag and drop the "Grip" block to the desired position in the programming field to grip a part.

If the part was successful gripped, the commands configured in the program branch "Holding" are executed. If the grip was not successful, the commands configured in the program branch "No_Part" are executed. In the case of a timeout or error, the corresponding program branches are executed.

12.3 Release a part (Release)

To release the gripped part, drag the "Release" block to the desired position in the programming field using drag & drop. Depending on the result of the command, the function block branches to "Released", "Timeout" or "Error".

ENGLISH
TM

TM

ENGLISH
TM

ENGLISH
TM

14 Translation of the original declaration of conformity

In terms of the EC Machinery Directive 2006/42/EG, Annex II, Part B

Manufacturer:

Weiss Robotics GmbH & Co. KG
Karl-Heinrich-Kaeferle-Str. 8
71640 Ludwigsburg, Germany

Distributor:

Weiss Robotics GmbH & Co. KG
Karl-Heinrich-Kaeferle-Str. 8
71640 Ludwigsburg, Germany

We hereby declare that the following products

GRIPKIT CR PLUS S, GRIPKIT CR PLUS L

meet the applicable basic requirements of the Machinery Directive (2006/42/EC). The incomplete machine may not be put into operation until conformity of the machine into which the incomplete machine is to be installed with the provisions of the Machinery Directive (2006/42/EC) is confirmed.

Applied harmonized standards, especially:


EN ISO 12100-1	Safety of machines – Basic concepts, general principles for design – Part 1: Basic terminology, methodology
EN ISO 12100-2	Safety of machines – Basic concepts, general principles for design – Part 2: Technical principles

The manufacturer agrees to forward the special technical documents for the incomplete machine to state offices on demand. The special technical documents according to Annex VII, Part B, belonging to the incomplete machine have been created.

Person responsible for documentation: Dr.-Ing. K. Weiß, Tel.: +49(0)7141/94702-0

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The specification in this document are subject to change without prior notice for the purpose of product improvement. All trademarks are the property of their respective owners

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Weiss Robotics GmbH & Co. KG, Management